


PRECISION COMFORT, QUICK-CLEAN STANDARD, VALVE CORE AND EURO-STYLE SYRINGE



ENGLISH

Syringe Parts

- #0084 8-32 Barb & Washer (Pkg. 10)
- #2201 O-Ring (Pkg. 12)
- #2202 O-Ring (Pkg. 12)
- #2203 O-Ring (Pkg. 12)
- #2211 O-Ring (Pkg. 12)
- #2212 O-Ring (Pkg. 12)
- #2267 O-Ring (Pkg. 12)
- #2293 O-Ring (Pkg. 12)
- #3060 Syringe Tip (1 ea.)
- #3180 Syringe Tip Indexing (1 ea.)
- #3061 Syringe Tip (Pkg. 5)
- #3181 Syringe Tip Indexing (Pkg. 5)
- #3059 Syringe Tip (Pkg. 50)
- #3182 Syringe Tip Index (Pkg. 50)
- #3065 Tip Adapter Assembly
- #3085 Button Valve Assembly
- #3093 Button Assembly, QC Syringe
- #3637 Button Valve Assembly
- #3090 Valve Core (Pkg. 10)
- #3491 Handle, Quick-Clean Syringe
- #9601 Retaining Ring
- #9648 Pin, Quick-Clean Syringe Button
- #9651 Pin, Euro-Style Syringe Button
- #9676 Spring, Syringe Button
- #9677 Spring, Tip Adapter



CLEANING

Clean the external surfaces of the syringe using a solution of mild detergent and warm water. A soft bristled brush may be used to clean around the buttons, tip collar and syringe threads.

Do not use powdered cleansers, scouring pads or abrasive scrubbers any of which can damage the surface finish of the syringe.

DISINFECTING

The use of an EPA-registered hospital disinfectant is adequate. Follow disinfectant manufacturer's instructions for use when disinfecting syringe heads.

TIP STERILIZATION

Clean tips under running water for 30 seconds with a soft brush and place in an ultrasonic bath with an enzymatic cleaner to remove superficial debris prior to sterilization; any debris may reduce the effectiveness of the sterilization. Sterilization process is by steam autoclave. Parameters for wrapped are 132°C (270°F) for 15 minutes at temperature with a 30-minute drying time when using a gravity displacement autoclave. Parameters for unwrapped are 132°C (270°F) for 3 minutes at temperature when using a gravity displacement autoclave.

MAINTENANCE

Syringe parts that are subject to wear under normal use are provided in a syringe repair kit included with your new syringe. Repair kit includes all necessary tools to repair the syringe.

REPLACEMENT PARTS

Use only genuine DCI replacement parts which are designed specifically for your syringe to properly function.

LUBRICANT

Use only silicone-based lubricant, DCI PN 8032, to lubricate syringe button assemblies and all other O-rings in the syringe.



CAUTION! The use of other types of lubricants may cause the O-rings deterioration resulting in syringe malfunction or failure.

DCI 3-way syringes are designed for professionals in the dental industry to supply air and water to the oral cavity during dental procedures for rinsing purposes. It is operated by depressing the air, water or both buttons simultaneously to precisely deliver water, mist or dry air for rinsing, cooling or drying the preparation area.

INSTALLATION

Turn off the air and water supplies. Determine which of the syringe tubes leads to the air button and connect that tube to the air outlet on the dental unit or syringe manifold. Connect the other tube to the water outlet. For circulating syringes, the third tube is the return water, and is routed to the drain. Recommended operating pressures: Air 80 psi, Water 40 psi.

OPERATION

The syringe head is imprinted with symbols indicating button functions. The button on the left is water and one on the right is air. Press both buttons simultaneously for steady uniform water spray.

CHANGING THE TIP

A sterile syringe tip should be installed for each patient. To prevent contamination, sterile tips should be handled with a sterile gauze pad. To remove the syringe tip, press down on the large collar. When you feel a soft "click", the tip may be pulled straight out. Hold the collar down and insert the new tip. Be sure to press it all the way in, and then release the collar.



CAUTION! Syringe tip may eject if it is not properly installed. Test the installation by giving a firm pull on the tip to assure that it is inserted all the way, and that the locking mechanism is fully engaged.

ASEPSIS

DCI dental syringes are considered a noncritical item. Use disposable barrier protection or a hospital disinfectant to disinfect the syringe between each patient or when the barrier has been compromised. When visibly contaminated with blood or other potentially infectious materials, an EPA-registered hospital disinfectant with a tuberculocidal claim should be used.

Syringe tips are considered semicritical and must be sterilized after each use.

TROUBLESHOOTING

Problem	Solution
Leakage around button	- Replace button assembly
Leakage at the syringe tip	- If water leaks, replace water button assembly - If air leaks, replace air button assembly - Replace small O-ring located behind the Quick-Change Adapter
Constant water spray when air button is pressed	- Replace tip - Replace small O-ring located behind the Quick-Change adapter
Leakage at syringe handle	- Verify barbed fittings are tightened - Replace nylon washer

STORAGE

Store in a dry, dust free location and do not mix with contaminated instruments.

Syringe tips shall be kept in a sterilization wrap or pouch to maintain sterility. Refer to the instructions for use from the sterilization packaging manufacturer to determine how long sterility can be maintained.

DISPOSAL

Contact your local, authorized dealer for disposal of this device to ensure compliance with local environmental regulations.

PRODUCT SUPPORT

For additional product support, contact DCI Customer Service at 1(800)624-2793 or visit www.dcionline.com

www.dcionline.com
800-624-2793
PN 92543, Rev. H, 2021-11-05

Applicable Patents:
U.S. Pat. No. 5,125,835
5,199,871



DCI International, LLC
305 N. Springbrook Road
Newberg, OR 97132 USA